

UPCERA



TT ONE PRE-SHADED ZIRCONIA

Reliable Materials Expert

TT ONE Pre-shaded

The pre-shaded zirconia is available in 16 VITA classical shades and 2 bleach shades. It is easy to operate with good stability and efficiency.



Strength
1000 MPa


Optimum Indications

- Full contour posterior crown
- Posterior bridge (3 units, 4 units)


Shades



Sintering program



Apply for



Open CAD/CAM System,
Zirkonzahn CAD/CAM System



Physical Characteristics

Density after sintering	$\geq 6.0 \text{g/cm}^3$
CTE(25-500°C)	$(10.5 \pm 0.5) \times 10^{-6} \text{K}^{-1}$
Accelerated aging surface monoclinic phase content	<5%
Chemical solubility after sintering	<100 $\mu\text{g/cm}^2$
Radioactivity	<0.1Bq/g
Sintering temperature	1480°C

Chemical Composition

$\text{ZrO}_2 + \text{HfO}_2 + \text{Y}_2\text{O}_3$	>96.5%
Y_2O_3	5.8%-9.7%
Al_2O_3	<0.5%
Fe_2O_3	<0.5%
Er_2O_3	<2.0%
Others oxides	<0.5%

Follow us on: